

SLIT INTERDIGITAL MICRO MIXER

SIMM GROUP CLASS SIMM-V2, HPIMM, SIMHEX, SSIMM

03




SIMM group class

Principle


This group class of micro mixers is a classic amongst all IMM chemical micro processing products. It has been used by a large number of customers, is cited multiple times in literature, and is indeed one of our best sellers.

They combine the regular flow pattern created by multi-lamination with geometric focussing which speeds up liquid mixing.

Due to this double-step mixing, the slit mixers are amenable to wide variety of processes such as mixing, emulsification, single-phase and multiphase organic synthesis. Extensive knowledge on hydrodynamics, mixing performance and reaction engineering for diverse applications of these mixers has been documented worldwide.



Interdigital flow passes slit to create multi-lamellae



Interdigital mixing principle

SIMM-V2



Individual parts of the SIMM-V2 device

Principle

This version has all the benefits of mixing using multi-lamination and focusing only. Deliberately avoiding volume expansion, the inner volume could be decreased to only 8 μl , coming along with improved fluidic connections, e. g. to pumps and tube reactors, as it employs HPLC connectors. Compared to the connectors of the standard version SSIMM, the HPLC joint to steel tubing improves leak tightness and higher pressure operation can be achieved.

Technical Data

Order number	SIMM-V2
Mixing principles	Multi-lamination
Size (L x B x H)	30 x 40 x 30
Connectors (Inlet/Outlet)	1/16" / 1/16" HPLC
Standard mixing channels (μm)	45 x 200
Standard material	Body: 1.4571 Inlay: 1.4435
Options	Other materials like Hastelloy, Monell or Titan on request

Operating Conditions

Temperature ($^{\circ}\text{C}$)	-40 – 220
Pressure stability (bar)	100
Flowrate (l/h)	0.04 – 2.5
Residence time (ms)	14.4 – 720
Inner volume (μl)	8
Max Viscosity (mPas)	10000
Leakage Class	< L _{0.001}

HPIMM



Individual parts of the HPIMM

Principle

This micro mixer was optimized using a metal sealing for tightening the two parts of the housing. As a consequence, the limits of pressure and temperature during operation are much higher than for flat-seal tightened devices. The mixer also comprises expansion-free outlet channel geometry, i.e. renounces on jet mixing, but relies on multi-lamination and geometric focusing only.

Technical Data

Order number	HPIMM
Mixing principles	Multi-lamination
Size (L x B x H)	25 x 21 x 37
Connectors (Inlet/Outlet)	1/16" / 1/16" HPLC
Standard mixing channels (μm)	25 x 21 x 37
Standard material	Body: 1.4571 Inlay: 1.4401
Options	Other materials like Hastelloy, Monell or Titan on request


Operating Conditions

Temperature ($^{\circ}\text{C}$)	-40 – 500
Pressure stability (bar)	600
Flowrate (l/h)	0.04 – 2.5
Residence time (ms)	27 – 1350
Inner volume (μl)	15
Max Viscosity (mPas)	10000
Leakage Class	< L _{0.001}

SLIT INTERDIGITAL MICRO MIXER


SIMM GROUP CLASS SIMM-V2, HPIMM, SIMHEX, SSIMM

SIMHEX



Slit Interdigital Mixer Heat Exchanger (SIMHEX)

SSIMM



Standard Slit Interdigital Micro Mixer (SSIMM)

Principle

This micro mixer was optimized considering a heat exchange function within the mixer, using a graphite sealing for tightening the two parts of the housing. As a consequence, the limits of pressure and temperature during operation are limited but conveniently provide the possibility of heating or cooling the device. The mixer also comprises expansion-free outlet channel geometry, i.e. renounces on jet mixing, but relies on multi-lamination and geometric focusing only.

Technical Data

Order number	SIMHEX
Mixing principles	Multi-lamination
Size (L x B x H)	25 x 25 x 20
Connectors (Inlet/Outlet)	1/16" / 1/16" HPLC
Standard mixing channels (µm)	40 x 300
Standard material	Body: 1.4571 Inlay: 1.4401
Options	Other materials like Hastelloy, Monell or Titan on request; incl. heat exchanger function

Operating Conditions

Temperature (°C)	-40 – 500
Pressure stability (bar)	50
Flowrate (l/h)	0.04 – 2.0
Residence time (ms)	18 – 900
Inner volume (µl)	10
Max Viscosity (mPas)	10000
Leakage Class	< L _{0,001}

Principle

This micro mixer is the classic one amongst all IMM chemical micro processing products. It combines the regular flow pattern created by multi-lamination with geometric focussing and subsequent volume expansion, which speeds up liquid mixing of the multi-lamellae and leads to jet mixing. Due to the volume expansion the mixer contains an inner volume of 40 µl and is only offered with non-stainless soft tube connectors.


Technical Data

Order number	SSIMM
Mixing principles	Multi-lamination
Size (L x B x H)	19 x 30 x 16.5
Connectors (Inlet/Outlet)	1/16" / 1/16" soft tube
Standard mixing channels (µm)	45 x 200
Standard material	Body: 1.4571 Inlay: 1.4435
Options	Other materials like Hastelloy, Monell or Titan on request

Operating Conditions

Temperature (°C)	-40 – 220
Pressure stability (bar)	3
Flowrate (l/h)	0.04 – 1.5
Residence time (ms)	72 – 3600
Inner volume (µl)	40
Max Viscosity (mPas)	10000
Leakage Class	< L _{0,001}

Slit mixer inlays SMI (for SIMM-V2 and SSIMM)



Laser ablation Inlay for SIMM-V2 and SSIMM


This inlay fits the standard mixer as well as the version 2.

For both versions the inlays have a size of 11.0 mm x 7.5 mm and ~ 3.6 mm thickness with different possible channel sizes and depths.

The following inlays are available:

- Laser-ablation (channel width 45 µm, 200 µm channel depth) made of stainless steel (SS 316L) as standard but other materials like Hastelloy etc. or other channel dimensions on request, order number SMI-Lasab45200
- LIGA technology (channel width 25 µm or 40 µm) made from silver or nickel on copper with 300 µm channel depth, order numbers SMI-Ni25, SMI-Ni40, SMI-Ag25 or SMI-Ag40
- ASE (thermally oxidised silicon, channel width 30 µm or 50 µm) with 100 µm channel depth, order numbers SMI-Si30 or SMI-Si50. As these inlays are only 0.6 mm thick, extra bases of 3.0 mm thickness are needed

SIMHEX inlays SMHXI



Laser-cut Inlay for SIMHEX


This inlay fits the slit interdigital mixer heat exchange exclusively.

The size of SMHXI inlays: 20 mm x 6 mm

The following inlays are available:

- Laser-cutted inlays (channel width 45 µm, 250 µm channel depth) made of stainless steel (SS 316L) as standard but other materials like Hastelloy etc. or other channel dimensions on request, order number SMHXI-45250

High pressure mixer inlays HPMI



Laser-cut Inlay for HPIMM

This inlay fits the high-pressure slit mixer exclusively.

The size of HPMI inlays: 8.0 mm in diameter and 250 µm in thickness

The following inlays are available:

- Laser-cutted inlays (channel width 45 µm, 250 µm channel depth) made of stainless steel (SS 316L) as standard but other materials like Hastelloy etc. or other channel dimensions on request, order number HPMI-Las45250

